


# Unlock Tissue Stiffness with Precision

*Advanced Shear Wave Elastography on the Vevo F2 Ultrasound Platform*



## Shear Wave Elastography (SWE)

on the Vevo F2 Imaging system offers unparalleled sensitivity for assessing tissue stiffness in small animal models.



## Scan the QR code

to learn more about shear wave elastography on our website.

## KEY FEATURES

### Fully Integrated on Vevo F2

Seamless integration for enhanced workflow allowing users to select 4 or 2 push pulses and easily manipulate data display.

### Quantitative Analysis & Quality Display Map

Measure shear wave velocity and stiffness in kPa for precise assessment; a quality display map provides insight into data acquisition quality.

### Real-time Visualization

Superimposed, real-time data display maps onto anatomical images for comprehensive analysis.

### Optimized for Small Animal Models

Tailored for high-resolution imaging of small animal tissues.

