

CYBERSECURITY ALERT

Dear FUJIFILM VisualSonics customer,

FUJIFILM VisualSonics, Inc. ("FF-VSI") is aware of the worldwide attack of the WannaCry ransomware worm and has diligently evaluated the potential impact on its ultrasound systems. After a thorough investigation, FF-VSI has determined that the Windows operating system underlying the VevoMD clinical imaging system contains a factory configured firewall that protects the system from the WannaCry ransomware worm.

FF-VSI is committed to providing the best in class services to its customers. Regardless of this exploit's ineffectiveness on the Vevo MD imaging system, to ensure that FF-VSI's customers can continue to rely on its products, FF-VSI has decided to develop a patch for the Microsoft vulnerability.

FF-VSI has additionally determined that the following pre-clinical products could be affected by the WannaCry ransomware worm and is in the process of developing a patch to correct the Microsoft operating system vulnerability. FF-VSI will notify customers when the patch is available.

Affected systems:

Vevo LAZR-X	Vevo 2100
Vevo LAZR	Vevo 3100
Vevo 1100	Vevo 770

Recommended preventive action for FF-VSI pre-clinical products:

Please take this opportunity to protect your saved data by backing it up to an offline location and remove your imaging system from the network until the provided patch can be installed.

Should you have any questions, please contact FF-VSI's technical support team:

Phone (U.S. or Canada) +1-877-657-8118
+1-888-482-9449 (toll free)

Phone (Outside U.S. or Canada)
+1-425-951-1330, or call your local representative

Fax +1-425-951-6700

E-mail support@visualsonics.com

Web www.visualsonics.com

Europe Service Center

Main: +31 20 751 2020

English support: +44 14 6234 1151

French support: +33 1 8288 0702

German support: +49 69 8088 4030

Italian support: +39 02 9475 3655

Spanish support: +34 91 123 8451

Asia Service Center +65 6380-5581